

Specification

Liquid version of the agar, recommended for studies of antibiotic MIC (Minimum Inhibitory Concentration).

Presentation

	Packaging Details	Shelf Life	Storage
20 Tubes Tube 16 x 113 mm with: 10 ± 0.3 ml	16x113 mm glass tubes, ink labelled, metal-Non injectable cap. - 20 tubes per box	12 months	2-25 °C

Composition

Composition (g/l):

Beef dehydrated infusion.....	2.00
Casein Hydrolysate Peptone.....	17.50
Starch.....	1.50
Ca++.....	0.02
Mg++.....	0.01

Description /Technique

Mueller Hinton Broth may be used for inoculum preparation for MIC tests and for preparation of antimicrobial dilutions for the microdilution or macrodilution procedure.

Note: For *S. pneumoniae*, use the medium with Lysed Horse Blood.

Quality control

Physical/Chemical control

Color : Pale yellow pH: 7.3 ± 0.2 at 25°C

Microbiological control

≤100 CFU. Min. 50 CFU (Productivity)

Analytical methodology according to ISO 11133:2014/A1:2018; A2:2020.

Aerobiosis. Incubation at 37 ± 1 °C, reading after 21 ± 3 h

Microorganism

Escherichia coli ATCC® 25922, WDCM 00013
Enterococcus faecalis ATCC® 29212, WDCM 00087
Ps. aeruginosa ATCC® 27853, WDCM 00025
Stph. aureus ATCC® 25923, WDCM 00034
Escherichia coli ATCC® 35218

Growth

Good
Good
Good
Good
Good

Sterility Control

Incubation 48 h at 30-35 °C and 48 h at 20-25 °C: NO GROWTH.
Check at 7 days after incubation in same conditions.

Bibliography

- BAUER, A.L., W.M.M. KIRBY, J.C. SHERRIS & M. TURCK (1966) Antibiotic susceptibility testing by a standardized single disc method. A. J. Clin. Pathol 45:493.
- BARRY, A.L., M.D. COYLE, C. THORNBERRY, E.H. GARLACH & R.W. HAWKINSON (1979) Methods of measuring zones of inhibition with Bauer-Kirby disk-susceptibility test. J. Clin. Microbiol. 10:885-889.
- CFR (1972) Rules and Regulations. 37:20525. Washington. DC. USA.
- CLSI (2006) Document M6-A2. Protocols for evaluating dehydrated Mueller-Hinton Agar: Approved Standard. 2nd ed. Clinical and Laboratory Standards Institute. Pennsylvania. USA.
- CLSI (2006) Document M2-A9. Performance standards for antimicrobial disk susceptibility tests: Approved Standard. 9th ed. Clinical and Laboratory Standards Institute. Pennsylvania. USA.
- ERICSSON & SHERRIS (1971) Acta Pathol. Microbiol. Scand Suppl 217 p:90.
- HINDLER, J. (1998) Antimicrobial Susceptibility Testing in Essential Procedures for Clinical Microbiology. ASM Press. Washington. DC. USA.
- MUNRO, S., (1995) Disk Difusión Susceptibility Testing, in Clinical Microbiology Procedures Handbook, H.D. Isenberg (ed) APHA. Washington. DC. USA.
- MILLER, J.M, C. THORNBERRY & C.N. BAKER (1984) Disk difusión susceptibility test troubleshooting guide. Lab. Med. 15:183-185.
- NEUMAN, M.A., D.F. SAMM, C. THORNSBERRY & I.E. MCGOWAN (1991) New developments in antimicrobial agent susceptibility testing: A practical guide. ASM. Washington. DC. USA.